THE BLACK


CROSS-LINKED HYALURONIC ACID
FILLER WITH LIDOCAINE



CE Approved Dermal Filler Formulated with High-Purity Hyaluronic Acid







THE BLACK

FINE

• HA: 20mg/mL

Lidocaine: 3mg/mL

· Volume: 1.0mL

Needle: 27G ½' / 30G ½

· Purpose: Fine Wrinkle Correction'

· Injection Depth: Mid Dermis



THE BLACK

DEEP

• HA: 20mg/mL

· Lidocaine: 3mg/mL

· Volume: 1.0mL

Needle: 27G ½' (2ea)

• Purpose: Deep Wrinkle Correction

· Injection Depth: Deep Dermis



THE BLACK

VOLUME

HA: 20mg/mL

· Lidocaine: 3mg/mL

· Volume: 1.0mL

• Needle: 25G ½' (2ea)

· Purpose: Volumizing & contouring

Injection Depth: Sub-Q Layer

Strong Points of THE BLACK HA Filler

Safety

Production under strict quality testing system dismisses traces of BDDE remain in the final product

Aesptic Manufacturing Environment

ISO13485, GMP, GMP facility certified

Smooth Injection

Stable extrusion force by fine-tuned gel particles

Superior Quality

Formulated with high-purity hyaluronic acid

Patients' Satisfaction

Unique HA structure provides harmonized look with long duration

Various Application

Design in various solutions, considering the characteristics of each treatment area







PLUMP



CONTOUR







MADE IN KOREA

H. www.wizmedi.kr Contact. wizmedisales@wizmedi.kr